

DrägerService®

Dräger

PMS Parts Reference Guide

Part Number: SP00216

Rev: -

Date: 31 May 2001

© 2001 Draeger Medical, Inc.

Microvent/Oxylog 2000

[RETURN TO SERVICE PROCEDURE TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

PMS Parts Reference Guide


This document is intended as a reference guide to assist in the replacement of routine maintenance parts.

Microvent / Oxylog 2000

<u>Description</u>	<u>Part Number</u>	<u>Comments</u> Microvent/Oxylog
12 month	4114768	
	R26807	O-Ring - 1x / 3x
	8410653	Diaphragm - 1x / 1x
	8400827	O-Ring - 2x / none
24 month	4114769	
	2M01524	Packing Ring - 1x / 2x
	D15330	O-Ring - 2x / 2x
	M16143	Washer - none / 2x
	8411599	Battery Pack - 1x / 1x
	D02316	Filter - 1x / 1x
	R24718	O-Ring - none / 1x
	2M12839	O-Ring - 1x / 1x
	R18566	O-Ring - 2x / none
72 month	<i>No kit</i>	
	8402745	Regulator - 1x / 1x
	<i>Or</i> 8411142	Reg Repair Set (Lorch)
	<i>Or</i> 8406678	Reg Repair Set (Knocks)


PMS Parts Reference Guide (continued)

Microvent Expiration Valve




PMS Parts Reference Guide (continued)

Oxylog 2000 Respiration Valve




PMS Parts Reference Guide (continued)

Microvent / Oxylog 2000 Flow Sensor




PMS Parts Reference Guide (continued)

Microvent




PMS Parts Reference Guide (continued)

Oxylog 2000



PMS Parts Reference Guide (continued)

Microvent / Oxylog 2000



[RETURN TO SERVICE PROCEDURE TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

[RETURN TO SERVICE PROCEDURE TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

DrägerService®

Dräger

DrägerService is a division of
Draeger Medical, Inc.
3122 Commerce Drive
Telford, PA 18969
Tel: (215) 721-5402
(800) 543-5047
Fax: (215) 721-5784
Web: www.draegermedical.com
Printed in the U.S.A.